

Molly[®] RNS

The Ideal Safety Solution



Technical Specifications

Device Type	Disposable Autoinjector
Primary Container	1 mL Long Pre-filled Syringe with RNS
Dosage Type	Fixed Full Contents
Needle Extension	Subcutaneous
Needle Insertion	Manual
Needle Removal	Manual
Activation Method	Needle cover

Handling Instructions



1.

Remove cap



2.

To activate, press
needle cover against
injection site & hold



+

Remove device
from injection site
& discard in sharps
container

Key Features

- ✓ Ultra compact design
- ✓ Needle cover activation
- ✓ Pre-configured or customized to client specifications
- ✓ Designed for Rigid Needle Shields